



**Cyrusbioscience**

## Total Nucleic Acid Isolation Kit-Virus

**Cat. Number:** TDRNAV

**Storage:** Store at Room Temperature.

**Sample:** Up to 200ul of serum, plasma, body fluids and the supernatant of viral infected cell cultures

### Kit Contents

| Total Nucleic Acid Isolation Kit-Virus | 100 Reactions |
|----------------------------------------|---------------|
| VN Column with 2 mL Collection Tube    | 100pcs        |
| Resuspension Buffer                    | 45mL          |
| Lysis Buffer (add Ethanol)             | 6mL (45mL)    |
| Wash Buffer 1                          | 45mL          |
| Wash Buffer 2 (add Ethanol)            | 15mL (60mL)   |
| Elution Buffer                         | 10mL          |
| This Datasheet                         | 1 copy        |

### Note:

- Add ethanol (96-100%) to Lysis Buffer and Wash Buffer 2, shake before use (see bottle label for volume).
- Check Buffers before use for salt precipitation. Re-dissolve any precipitate by warming to 37°C.
- Resuspension Buffer and Wash Buffer 1 contain irritants. Wear gloves when handling these buffers.

### Additional requirements

- Absolute ethanol
- PBS (Phosphate Buffered Saline)
- Microcentrifuge tubes (DNase and RNase free)

### Description


The Total Nucleic Acid Isolation Kit-Virus provides a fast, simple, and cost-effective method for isolation of viral DNA/RNA from cell-free samples such as serum, plasma, body fluids and the supernatant of viral infected cell cultures. Its unique buffer system efficiently lyses cells and degrades proteins, allowing for nucleic acids to bind to the glass fiber matrix of the columns easily. Contaminants such as salts, metabolites and soluble macromolecular cellular components are removed through the Wash step. Phenol extraction and ethanol precipitation are not required, and high-quality Nucleic Acids are eluted with RNase-free elution buffer. Viral DNA/RNA isolated with Total Nucleic Acid Isolation Kit-Virus is suitable for a variety of routine applications, including Real-time PCR/RT-PCR, Automated Fluorescent DNA Sequencing, PCR, and other enzymatic reactions. The entire procedure can be completed within 15-20 minutes.



**Cyrusbioscience**

---

## Procedure for Preparation of Total Nucleic Acid Isolation



### Step 1. Lysis

1. Transfer up to 200  $\mu$ l of virus sample into a 1.5 ml microcentrifuge tube and add 400  $\mu$ l of Resuspension Buffer. (If the sample is less than 200  $\mu$ l, adjust the sample volume to 200  $\mu$ l with PBS).
2. Mix well and let stand at room temperature for 10 minutes.  
# Pre-heat the Elution Buffer to 75°C for Step 4 DNA Elution.

### Step 2. Nucleic Acid Binding

1. Add 450  $\mu$ l of Lysis Buffer (ethanol added) to the sample lysate and shake vigorously.
2. Place a VN Column in a Collection Tube. Transfer 700  $\mu$ l of the lysate mixture to the VN Column.
3. Centrifuge at 16,000 x g for 1 minute. Discard the flow-through and place the VN Column back in the Collection Tube.
4. Transfer the remaining lysate mixture to the VN Column.
5. Centrifuge at 16,000 x g for 1 minute. Discard the flow-through and place the VN Column back in the Collection Tube.

### Step 3. Wash

1. Add 400  $\mu$ l of Wash Buffer 1 into the VN Column. Centrifuge at 16,000 x g for 30 seconds. Discard the flow-through and place the VN column back into the Collection tube.
2. Add 600  $\mu$ l of Wash Buffer 2 (ethanol added) into the VN Column. Centrifuge at 16,000 x g for 30 seconds. Discard the flow-through and place the VN column back in the Collection tube.
3. Centrifuge at 16,000 x g again for 2 minutes to remove residual Wash Buffer 2 Buffer.

### Step 4. Elution

1. Place the VN column in a clean 1.5 ml microcentrifuge tube (DNase and RNase free).
2. Add 50-100  $\mu$ l of Pre-heated Elution Buffer or RNase-free water (pH is between 7.0 and 8.5) to the center of each VN column, let stand for 2 min, and centrifuge at 14,000 x g for 2 min.



**Cyrusbioscience**

## Troubleshooting

| Problems                                           | Possible Causes                                                        | Suggestions                                                                                                                                                                     |
|----------------------------------------------------|------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Poor performance of RNA in downstream applications | Interference of the residual ethanol                                   | Be sure to remove Lysis Buffer and Wash Buffer 2 completely.                                                                                                                    |
| Low yields                                         | Insufficient performance of the elution buffer during the elution step | Remove the residual buffers during the wash steps completely. These residual buffers decrease the efficiency of the following elution steps.                                    |
|                                                    | Incomplete lysis                                                       | Check the incubation time of the Lysis Step.                                                                                                                                    |
|                                                    | Viral nucleic acid remains on the column                               | Repeat the Elution Step with the eluant.<br>Extend the time Elution Buffer or RNase-free water stays on column centers from Step 4 from 2 to 5 minutes prior to centrifugation. |
| Degraded RNA                                       | Source                                                                 | Do not freeze and thaw sample more than once.<br>Increase the virus concentration in the sample.                                                                                |
|                                                    | RNase contamination                                                    | Be sure not to introduce RNase during the procedure.<br>Check buffers for the RNase contamination.                                                                              |