



Cyrusbioscience

GeneKlean® Total RNA Isolation Kit-Plant

Cat. Number: TRNAP

Storage: Store at Room Temperature.

Sample: 100mg of fresh plant tissue or 25mg of dry plant tissue

Yield: Up to 30ug Total RNA

Kit Contents

Total RNA Isolation Kit-Plant	100 Reactions
GeneKlean® Total RNA Column with 2 mL Collection Tube	100
Lysis Buffer	110mL
Wash Buffer 1	45mL
Wash Buffer 2 (add Ethanol)	15mL (60mL)
Elution Buffer	10mL
This Datasheet	1 copy

Note

- Always wear latex or vinyl gloves while handling reagents and RNA samples to prevent RNase contamination.
- Add ethanol (>99.5%) to Wash Buffer 2, shake before use (see bottle label for volume).
- Check Buffers before use for salt precipitation. Redissolve any precipitate by warming to 37°C.
- Buffers contain irritants. Wear gloves when handling these buffers.

Additional requirements

- Ethanol (>99.5%)
- Isopropanol
- RNase-free pipet tips and 1.5 ml microcentrifuge tubes
- b-mercaptoethanol

Optional Step for DNA Degradation:

Add 2 µl DNase I (2KU/ml) and 10 µL reaction buffer (300 mM Tris-HCl (pH 7.5), 60 mM MnCl₂, 300 µg/ml BSA) to the 50µl final product. Let stand for 10 minutes at room temperature (at 25°C).


Description



Cyrusbioscience

GeneKlean® Total RNA Isolation Kit-Plant provides a fast, simple, and cost-effective method for isolation of total RNA from plant samples. Detergents and chaotropic salt are used to lyse cells and inactivate RNase. The specialized high-salt buffering system allows RNA species longer than 100 bases to bind to the glass fiber matrix of the spin column. GeneKlean® Total RNA Isolation Kit-Plant is suitable for a variety of routine applications including RT-PCR, cDNA Synthesis, Northern Blotting, Differential display, Primer Extension and mRNA Selection. The entire procedure can be completed within 60 minutes.

Procedure for Preparation of Plant Total RNA Isolation



Sample Preparation

1. Cut off 100 mg of fresh plant tissue or 50 mg of dry plant tissue.
2. Grind the sample under liquid nitrogen to a fine powder using a mortar and pestle.

Step 1. Lysis

1. Add 1mL of Lysis Buffer and 10uL of b-Mercaptoethanol to the sample in the motor and grind the sample until it is dissolved completely.
2. Transfer the dissolved sample to RNase-free 1.5mL microtube. Incubate at 70°C for 30 mins.(Invert the tube every 10 mins)
3. Centrifuge at 4°C at 15,000 x g for 10 mins. Transfer the supernatant to a new 1.5mL microtube.

Step 2. RNA Binding

1. Add 1/2 volume of Isopropanol to the sample from Step 1 and shake vigorously (eg. 250uL of Isopropanol to 500uL of sample)
2. Place GeneKlean® Total RNA Column in a 2mL Collection Tube.
3. Transfer the sample mixture to the GeneKlean® Total RNA Column. Centrifuge at 15,000x g for 30 secs.
4. Discard the flow-through and transfer the remaining mixture to the same GeneKlean® Total RNA Column.
5. Centrifuge at 15,000x g for 30 secs.
6. Discard the flow-through and place the GeneKlean® Total RNA Column back to the 2mL Collection Tube.

Step 3. Washing



Cyrusbioscience

1. Add 400uL of Wash Buffer 1 to the GeneKlean® Total RNA Column. Centrifuge at 15,000x g for 30 secs.
2. Discard the flow-through and place the GeneKlean® Total RNA Column back to the 2mL Collection Tube.
3. Add 600uL of Wash Buffer 2(ethanol added) into GeneKlean® Total RNA Column. Centrifuge at 15,000x g for 30 secs.
4. Discard the flow-through and place the GeneKlean® Total RNA Column back to the 2mL Collection Tube.
5. Centrifuge again for 3 mins at 15,000x g to dry the column matrix.

Step 4. Elution

1. Place GeneKlean® Total RNA Column. Centrifuge in a clean RNase Free 1.5mL microtube.
2. Add 50-200uL of Elution Buffer to the center of GeneKlean® Total RNA Column, stand for 2 mins. Centrifuge 2 mins at 15,000x g

Troubleshooting

Problems	Possible Causes	Suggestions
Degraded RNA/Low integrity	RNase Contamination	Clean everything, use barrier tips, wear gloves and a lab coat, and use RNase-free enzymes, Ex: RNase inhibitor.
Low yields of RNA	Incomplete lysis and homogenization	Don't use more samples than the suggested limit.
	Incorrect elution conditions	Add 100 uL of the Elution Buffer to the center of each Column, let it stand for 2 minutes, and centrifuge at 15,000 x g for 2 minutes.
Inhibition of downstream enzymatic reactions	Presence of ethanol in the purified RNA	Repeat the wash step: Centrifuge at 15,000 x g again for 2 minutes to remove the residual of Wash Buffer 2.