




Cyrusbioscience

Specification

DMSO, Dimethyl Sulfoxide



Product Code:	101-67-68-5
CAS number:	67-68-5
Formula:	(CH ₃) ₂ SO
Formula Weight:	78.13
Storage:	Store at RT

Specification

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Clarity)	Clear
Assay	>99.5%
Identity (IR)	Pass
pH	6~8
Purity (G.C.%)	≥ 99.9%
UV Spectroscopy	Pass
Water Content (K.F.)	≤ 0.1%
Reducing Material (As Formaldehyde)	≤ 0.01%

Responsible Laboratory Manager Quality Control: Taylor Lin

This document has been produced electronically and is valid without a signature.