




Cyrusbioscience

Specification

EDTA, Disodium Salt, Dihydrate

Ethylenediaminetetraacetic Acid, 2Na, 2H₂O



Product Code:	101-6381-92-6
CAS number:	6381-92-6
Formula:	C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ ·2H ₂ O
Formula Weight:	372.24
Storage:	Store at Room Temperature. Keep Dry.

Specification

Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
160 mg/mL, 3 M NaOH	
Purity, %	>=99.0
Chloride(Cl ⁻), %	<=0.01
Sulfate(SO ₄ ²⁻), %	<=0.05
Iron(Fe), %	<=0.001
Heavy metal(as Pb), %	<=0.001

Responsible Laboratory Manager Quality Control: Taylor Lin

This document has been produced electronically and is valid without a signature.