




Cyrusbioscience

---

## Specification

### Tris HCl

(Tris-(Hydroxymethyl)-aminomethane hydrochloride)



Product Code:	101-1185-53-1
CAS number:	1185-53-1
Formula:	C <sub>4</sub> H <sub>12</sub> NO <sub>3</sub> Cl
Formula Weight:	157.60
Storage:	Store at RT.

### Specification

Appearance:	White Crystalline Granule
Assay, %:	>=99.0%
Moisture(LOD):	<=1.0%
pH (0.5M in H <sub>2</sub> O):	3.5-5.0
UV Abs., 260nm, 0.5M in H <sub>2</sub> O	<=0.015
UV Abs., 280nm, 0.5M in H <sub>2</sub> O	<=0.010

Responsible Laboratory Manager Quality Control: Taylor Lin

This document has been produced electronically and is valid without a signature.