




Cyrusbioscience

Specification

X-GlcA

5-Bromo-4-chloro-3-indolyl β -D-glucuronide sodium salt



Product Code:	101-129541-41-9
CAS number:	129541-41-9
Formula:	C ₁₄ H ₁₂ BrClNO ₇ Na
Molecular mass.:	444.6
Storage:	Store at 4°C . Protect from light and humidity. For long term storage hold at -20°C .

Specification

Appearance:	White powder
Purity(TLC):	>98%
Identity(¹ H-NMR):	Conforms to structure

Responsible Laboratory Manager Quality Control: Taylor Lin

This document has been produced electronically and is valid without a signature.