



Cyrusbioscience

## Cy2Ab® Goat Anti Mouse IgG(H+L)-HRP

### Secondary Antibody

**Product number:** 2AB0001

**Package:** 500ul

**Description:** Goat anti-Mouse IgG (H+L)-HRP is conjugated with high-activity HRP.

**Product type:** Polyclonal

**Concentration:** Lot-dependent, please refer to product label.

**Immunogen:** Purified mouse IgG

**Host:** Goat

**Reactivity:** Mouse (Minor cross-reactivity to rat)

**Conjugation:** HRP (RZ>3)

**Isotype:** N/A

**Purity:** Immunogen affinity purified

**Form:** Liquid

**Storage buffer:** Phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, and 50% glycerol, pH 7.4.

**Storage:** Store at 4°C for at least 6 months.

Store at -20°C for long term preservation.

Avoid repeated freeze / thaw cycles.

**Suggestion:** Divide antibody into several vials. Keep only using vials at 4°C.

**Application:** WB, ELISA, IHC

**Manual:** 1:5,000-1:25,000 for WB and ELISA; 1:500-1:2,500 for IHC

Working dilution for specific application should be determined to obtain the best conditions.